

# BRUNO

# Wheelchair-Lifter



- ◆ The Bruno Wheelchair-Lifter Model AWL-100 Vehicle Lift eliminates the struggle of lifting/storing your folding wheelchair into the trunk of your car or station wagon. You don't need to struggle any longer to get your wheelchair in/out of the trunk of your car.
- ◆ The Wheelchair-Lifter Model AWL-100 raises and tilts your wheelchair horizontally, all at the touch of a button. Plus, Bruno has two additional Wheelchair-Lifter Models – the AWL-150 rear-door application and the AWL-160 side-door application – to vertically lift and store your wheelchair into a full-size van or minivan.
- ◆ This electrically-operated lift is powered by the vehicle's battery; it can also be used with Bruno's optional Pow'rPac . When the wheelchair is above the trunk, simply swing the lift arm and wheelchair in and push the **DOWN** button; then fold the lift arm down.
- ◆ The Wheelchair-Lifter Model AWL-100 can be installed in most auto trunks, station wagons, vans, minivans, sport utility vehicles (SUVs) and hatchbacks as long as the storage area is large enough to accommodate your folding wheelchair. Please see the back of this brochure for a step-by-step illustration of how the Wheelchair-Lifter Model AWL-100 works.

## OPTIONAL BRUNO POW'RPAC



- ◆ Easy to use!
- ◆ Convenient!
- ◆ Works in most minivans, SUVs, station wagons & sedans!

The complete, independent power supply for your automotive lift!

Pow'rPac easily installs in your vehicle – just plug into the vehicle's power outlet!

- ◆ **EXCLUSIVE!** Bruno's R<sup>3</sup> Lift Belt is 50% stronger! Count on Bruno's R<sup>3</sup> Lift Belt for ease of use, longer life and safer operation – *another Industry First!*
- ◆ Expert installation and complete, after-sale service is provided by Bruno's worldwide network of dealers.
- ◆ Only Bruno has over 25 different interior/exterior scooter, wheelchair and power chair lifts for nearly any sedan, minivan, sport utility vehicle, full-size van, pickup truck and the new cross-over vehicles.
- ◆ All Bruno lifts are eligible for reimbursement under the terms of Mobility Programs offered by GM, Saturn, DaimlerChrysler/Dodge, Ford, Toyota, Lexus, Volvo and Volkswagen.



# Using the *Bruno Wheelchair-Lifter*<sup>®</sup> is easy!



Once the wheelchair is positioned at the rear of your car, remove the seat and hook the *lift claw* to the T-bar docking device.



To lift, simply push the **UP** button. The lift will tilt the wheelchair into the horizontal position with your guidance.



Continue to push the **UP** button until the wheelchair is above the trunk opening.



Swing the wheelchair above the trunk opening.



Push the **DOWN** button to lower the wheelchair into the trunk and remove the *lift claw* from the docking device. The lift and wheelchair fit snugly in your trunk.



**WHEELCHAIR-LIFTER Model AWL-150** REAR-DOOR APPLICATION  
**WHEELCHAIR-LIFTER Model AWL-160** SIDE-DOOR APPLICATION



The Lift weight rating is the **MAXIMUM** load value. Factors inherent to specific applications including, but not limited to, vehicle structure and condition, operating environment, and customer ability may reduce the effective weight rating.

**Bruno lifts are eligible for reimbursement under the terms of Mobility Programs offered by GM, Saturn, DaimlerChrysler/Dodge, Ford, Toyota, Lexus Volvo and Volkswagen.**

**Represented By:**



Bruno Independent Living Aids, Inc., 1780 Executive Drive,  
 P.O. Box 84, Oconomowoc, WI 53066  
 1-262-567-4990 ♦ FAX: 1-262-953-5501  
[www.bruno.com](http://www.bruno.com)